

Microcide TB TECHNICAL BULLETIN

Cleaner
Deodorizer
Disinfectant
Virucide
Tuberculocidal
Fungicide
Mildewstat
Fast Kill Times



Microcide TB Kills 99.9%* of SARS-CoV-2, which causes COVID-19 on hard non-porous surfaces in just 1 minute!

Proven Effective Against:

CA-MRSA,
MRSA,
VRE,
VISA,
MRSE,
HIV-1,
HAV,
HBV & HCV,
TB,
SARS-COV-2,
Canine Parvovirus,
Norovirus,
(and many others, see inside for details)

Meets OSHA Bloodborne Pathogen
Standard for HIV, HBV & HCV.

EPA Registration Number
1839-83-5449

EPA Establishment Number
8155-OH-1

US Patent No. 5,444,094

SPECIFICATIONS FOR MICROCID TB

Appearance Clear Blue Liquid
Fragrance Citrus
Density lbs./US gal.(g/ml)....8.38 (1.01)
pH 11.7
Contact time:.....30 Seconds to 10 minutes
(see label for specific instructions)
Viscosity cps @ 25°C.....< 100
Flashpoint SETA °C (°F)...> 94 (>201)
RVOC, US EPA%0
Phosphate (%).....0
EPA Number 1839-83-5449
EPA Est. Number.....8155-OH-1

Microcide TB Disinfectant Cleaner

Microcide TB Disinfectant Cleaner is designed specifically as a non-acid, intermediate level, ready-to-use disinfectant and cleaner for use in hospitals, nursing homes, patient rooms, shower rooms, locker rooms, public restrooms, schools, office buildings, and athletic/recreational facilities.

Microcide TB is formulated to disinfect hard, non-porous, inanimate environmental surfaces: floors, walls, stainless steel surfaces, glazed porcelain, glazed ceramic tile, plastic surfaces, chrome, brass, copper, laminated surfaces, bathrooms, shower stalls, bathtubs, and cabinets. For plastic and painted surfaces, spot test on an inconspicuous area before use.

Microcide TB is an intermediate level disinfectant, meaning it has a claim against TB (*Mycobacterium bovis*). In addition to this claim, Microcide TB also has claims for several non-enveloped viruses.

Non-Enveloped Viruses:

Non-enveloped viruses differ from the more common enveloped viruses in their method of transmission between cells. These differences contribute to non-enveloped viruses being much more difficult to eliminate than enveloped viruses, in terms of cleaning and disinfecting. Standard disinfecting products and procedures sometimes do not efficiently disinfect against non-enveloped viruses.

Microcide TB meets the CDC criteria for disinfectant products with label claims for non-enveloped viruses: Norovirus, Adenovirus Type 2.

Microcide TB is intended for use on hard, non-porous surfaces, follow label instructions for non-enveloped viruses.

Microcide TB has claims against the following tough-to-eliminate non-enveloped viruses:

- **Canine Parvovirus** is a highly contagious virus that affects dogs and is therefore of high concern in veterinary facilities. Treatment often requires hospitalization, and mortality rates can reach up to 91% in untreated cases.
- **Norovirus** is a highly contagious virus that causes nausea, vomiting, and diarrhea. This is often called “stomach flu” or “the cruise ship virus.” Although severe illness is rare, this virus can run rampant when an infected person is in contact with many others, such as on cruise ships or in schools.
- **Rhinovirus** is one of the most common infectious agents among humans, a cause of the common cold. Symptoms vary and include nausea, headache, cough, sore throat, and fever.
- **Rotavirus** is particularly concerning among children, and can often lead to hospitalization due to dehydration. Symptoms include fever, vomiting, and diarrhea.

Dwell Times: Microcide TB has rapid kill times for a number of pathogens, including: MRSA, VRE, and VISA (3 minutes), Influenza (2 minutes), and Norovirus (30 seconds).

Microcide TB Disinfects hard non-porous surfaces by killing 99.9% of SARS-CoV-2 the Novel Coronavirus that causes COVID-19 in just One Step.

Multi-Clean Microcide TB Disinfectant Cleaner (EPA registration No. 1839-83-5449) has demonstrated effectiveness against viruses similar to the 2019 Novel Coronavirus (SARS-CoV-2, the cause of COVID-19) on hard, non-porous surfaces.

Therefore, **Multi-Clean Microcide TB Disinfectant Cleaner** can be used against the 2019 Novel Coronavirus (SARS-CoV-2, the cause of COVID-19) when used in accordance with directions for use against Canine Parvovirus on hard, non-porous surfaces.

MICRO-ORGANISMS

BACTERICIDAL ACTIVITY (BACTERIA): Using AOAC test methods under Good Laboratory Practices, in the presence of 5% blood serum and a 3 minute contact time, unless otherwise noted, Microcide TB kills the following organisms on hard non-porous inanimate surfaces:

Organism	Contact Time
Pseudomonas aeruginosa	3 minutes
Staphylococcus aureus	3 minutes
Salmonella (choleraesuis) enterica	3 minutes
Escherichia coli	3 minutes
Escherichia coli O157:H7	3 minutes
Streptococcus pyogenes (Necrotizing Fasciitis-Group A) (V.A. Medi-cal Center Isolate 04001)	3 minutes
Listeria monocytogenes	3 minutes
Yersinia enterocolitica	3 minutes
Enterococcus faecium	3 minutes
Corynebacterium ammoniagenes	3 minutes
Salmonella (typhi) enterica	3 minutes

ANTIBIOTIC-RESISTANT:

Using AOAC test methods under Good Laboratory Practices, in the presence of 5% blood serum and a 3 minute contact time, unless otherwise noted, Microcide TB kills the following organisms on hard non-porous inanimate surfaces:

Methicillin resistant Staphylococcus aureus (MRSA)	3 minutes
Vancomycin resistant Enterococcus faecalis (VRE)	3 minutes
Vancomycin intermediate resistant Staphylococcus aureus (VISA) (CDC Isolate 99287)	3 minutes
Methicillin resistant Staphylococcus epidermidis (MRSE)	3 minutes
Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA) (NRS 123) Geno-type USA400	3 minutes
Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA) (NRS 384) Geno-type USA300	3 minutes

VIRUCIDAL ACTIVITY (Viruses)

Contact Time

Virus Enveloped	
HIV-1 (associated with AIDS)	1 minute
Avian Influenza A strain H3N2	2 minutes
Avian Influenza A strain H9N2	2 minutes
Hepatitis B Virus (HBV)	1(or)5minute(s)
Hepatitis C Virus (HCV)	1(or)5minute(s)
Duck Hepatitis B Virus (DHBV)	1(or)5minute(s)
Bovine Viral Diarrhea Virus (BVDV)	1(or)5minute(s)
Human Coronavirus	2 minutes
SARS associated Coronavirus	2 minutes
SARS-Related Coronavirus 2 [SARS-CoV-2] USA-WA1/2020 BEI NR-52281	1 minute
Rabies Virus	30 seconds
Paramyxovirus (Mumps)	3 minutes
Virus Large Non-Enveloped	
Rotavirus	3 minutes
Virus Small Non-Enveloped	
Norovirus (Norwalk virus), Feline Calicivirus	30 seconds
Canine Parvovirus	10 minutes
Poliovirus Type 1	10 minutes
Rhinovirus type 39	3 minutes
Hepatitis A Virus (HAV) (University of Ottawa)	10 minutes

Mildew / Fungi

Contact Time

Trichophyton interdigitale, formerly Trichophyton mentagrophytes (athlete's foot fungus)	10 minutes
---	------------



600 Cardigan Road
Shoreview, Minnesota
Phone 651-481-1900