

# Microcide TB

## Disinfectant Cleaner



Microcide TB Disinfectant Cleaner is designed specifically as an intermediate level disinfectant for use in hospitals, nursing homes, patient rooms, shower rooms, locker rooms, public restrooms, schools, office buildings, and athletic/recreational facilities, where housekeeping is of prime importance in controlling the hazard of cross contamination. Microcide TB offers a faster kill time than many disinfectants for common pathogens.

**Tuberculocidal:** Microcide TB kills **TB (Mycobacterium tuberculosis)** in **5 minutes**.

**Hospital Grade:** Microcide TB kills **Pseudomonas aeruginosa, Staphylococcus aureus, and Salmonella enterica (hospital trio)** in **3 minutes**.

**Antibiotic-Resistant Bacteria:** Microcide TB kills **MRSA, CA-MRSA, VRE, VISA, and MRSE** in **3 minutes**.

**Non-Enveloped Viruses:** Microcide TB is effective against the following non-enveloped viruses:

- **Poliovirus (10 minutes)**
- **Rotavirus (3 minutes)**
- **Norovirus (30 seconds)**

**Norovirus:** Microcide TB is effective against **Norovirus**, the rampant “cruise ship” virus that causes stomach flu, in **30 seconds**.

**Other Viruses:** Microcide TB is also effective against many other important viruses such as **Hepatitis A, B and C, Influenza A, H1N1, and Canine Parvovirus**.

### Kills Norovirus in 30 seconds!

Microcide TB meets the CDC criteria for disinfectant products with label claims for non-enveloped viruses: Norovirus, Adenovirus Type 2.

Microcide TB is intended for use on hard, non-porous surfaces, follow label instructions for non-enveloped viruses.



Microcide TB can be used in the E-Spray Electrostatic Sprayers.



Meets OSHA Bloodborne Pathogen Standard for HIV, HBV, and HCV. US Patent No. 5,444,094

### Properties

Appearance .....	Light Blue/Green Liquid
Fragrance .....	Citrus Scent
pH .....	11.0 - 12.0
Activity .....	2,100 ppm Quat
Contact Time .....	30 Seconds to 10 Minutes (Read label for specific times needed to kill specific organisms.)
Dilution .....	Use Straight, no dilution
Flashpoint .....	>201 F
Storage Conditions .....	40 F - 110 F (4 C - 43 C)
Storage Stability .....	Excellent, 1 year minimum
Freeze-Thaw Stability .....	Will freeze but is not damaged. Agitate after thawing.
EPA Registration Number .....	1839-83-5449
VOC Compliant - VOC Free	

### Packaging

12 qt. case, 910248  
5 gallon, 910245  
55 gallon, 910243



600 Cardigan Road  
Shoreview, Minnesota  
Phone 651-481-1900  
Web [www.multi-clean.com](http://www.multi-clean.com)